



## Axolotl DualGraft™

Indicated as a wound covering and structural barrier

Processed through minimally manipulated techniques **designed to retain the qualities of the native ECM**

**Simplified application** | Processed with **epithelial surfaces facing outwards**

**NATURAL**  
**EFFICIENT**  
**EFFECTIVE**



## Evidence supporting the use of Amnion

### The Extracellular Matrix (ECM)

3-D structure & essential component in all human tissues that aids in: Cellular adherence | Tissue anchoring  
Cellular signaling | Recruitment of cells<sup>1</sup>

### ECM contributes to:

Formation of granulation tissue | Germination of new blood vessels | Maintaining viability of new tissue<sup>2</sup>

### Human Amniotic Membrane (hAM)

Surgical recommendations are for application of hAM epithelial side against the wound<sup>3,4,5</sup>




## Axolotl DualGraft™ Application



1  
Open box, remove outer tray. Open the outer tray and remove the inner tray (the inner tray is sterile and may be placed onto the sterile field). Open the inner tray and remove the Axolotl DualGraft package.



2  
Using sterile scissors, cut open the package and remove Axolotl DualGraft.



3  
Apply Axolotl DualGraft directly to the area of interest. Secure Axolotl DualGraft with sutures, glue, or steri-strips if desired, based on physician preference. Cover graft with a wide meshed non-adherent dressing (Adaptic). Use a bolster dressing such as gauze if necessary to fill dead space.



**STORAGE**  
Stored at **ambient temperature**



**PRESERVATION METHOD**  
Terminal irradiation gives the membrane a **5-year shelf life**.



**CONTAINS ESSENTIAL PROTEINS & GROWTH FACTORS**  
Proteins found in Axolotl DualGraft™ include:  
• Collagen I, III, IV, V and VII  
• Fibronectin

| SKU     | DESCRIPTION        | SIZE  | Q CODE |
|---------|--------------------|-------|--------|
| ADG12   | Axolotl DualGraft™ | 1x2   | Q4332  |
| ADG22   | Axolotl DualGraft™ | 2x2   | Q4332  |
| ADG24   | Axolotl DualGraft™ | 2x4   | Q4332  |
| ADG44   | Axolotl DualGraft™ | 4x4   | Q4332  |
| ADG46   | Axolotl DualGraft™ | 4x6   | Q4332  |
| ADG48   | Axolotl DualGraft™ | 4x8   | Q4332  |
| ADG1010 | Axolotl DualGraft™ | 10x10 | Q4332  |

Axolotl DualGraft has received a positive recommendation from the FDA Tissue Review Group (TRG) and is manufactured in an AATB accredited facility.